

Product datasheet for **MP208361**

Mrtfa Mouse qPCR Primer Pair (NM_153049)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	223701
Forward Sequence:	GCCTAGCCAATGGAAGCATCGT
Reverse Sequence:	CACCTTTGGCTTCAGCTCCTTG
ACCN:	BC050941 , NM_153049 , NM_153049.1 , NM_153049.2 , NM_153049.3 , BB850443 , BC054801 , BY200528 , BY208699 , BY749697 , BY750304
UniProt ID:	Q8K4J6
Synonyms:	AMKL; Bsa; Mal; Mrtf-A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)