

Product datasheet for **MP208402**

Mospd1 Mouse qPCR Primer Pair (NM_027409)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70380
Forward Sequence:	CAGTGAAGCCACAGTGCTGTGT
Reverse Sequence:	GGCTTTGTTTCAGAAACCTGGAGG
ACCN:	<u>BC018329</u> , <u>BC085100</u> , <u>BC100442</u> , <u>NM_027409</u> , <u>NM_027409.1</u> , <u>NM_027409.2</u> , <u>NM_027409.3</u> , <u>NM_027409.4</u> , <u>NM_027409.5</u> , <u>BF467665</u> , <u>BG068741</u> , <u>BG077592</u>
UniProt ID:	<u>Q8VEL0</u>
Synonyms:	1810018L05Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)