

## Product datasheet for **MP208449**

### Msi2 Mouse qPCR Primer Pair (NM\_054043)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	76626
Forward Sequence:	CAGCCGAAAGAAGTCATGTTCCC
Reverse Sequence:	CCTCTGCCATAGGTTGCCACAA
ACCN:	<u><a href="#">BC109368</a></u> , <u><a href="#">BC111809</a></u> , <u><a href="#">NM_054043</a></u> , <u><a href="#">NM_054043.1</a></u> , <u><a href="#">NM_054043.2</a></u> , <u><a href="#">NM_054043.3</a></u> , <u><a href="#">BK001483</a></u> , <u><a href="#">BQ174763</a></u>
UniProt ID:	<u><a href="#">Q920Q6</a></u>
Synonyms:	1700105C15Rik; AI451722; AI563628; AW489193; Msi2h
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)