

Product datasheet for **MP208454**

Mbnl2 Mouse qPCR Primer Pair (NM_175341)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	105559
Forward Sequence:	AGTCCCAACAGAGGTTCTACCC
Reverse Sequence:	GGCACAGTTTCCTCGCTGGAAC
ACCN:	<u>BC075665</u> , <u>NM_175341</u> , <u>NM_175341.1</u> , <u>NM_175341.2</u> , <u>NM_175341.3</u> , <u>NM_175341.4</u> , <u>BB673209</u> , <u>BC020406</u> , <u>BC023895</u> , <u>BC026711</u> , <u>BC079898</u> , <u>BE336239</u> , <u>BM932374</u> , <u>BQ946035</u> , <u>BY198796</u>
UniProt ID:	<u>Q8C181</u>
Synonyms:	1110002M11Rik; AI047808; AI837313; AI849185; AL118326; MBLL; mKIAA4072; R75232
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)