

Product datasheet for **MP208469**

Mybbp1a Mouse qPCR Primer Pair (NM_016776)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18432
Forward Sequence:	AACCCAGCAACTTGATGCGCCA
Reverse Sequence:	CTCATCCTCGTTTGTACAGGG
ACCN:	<u>BC048858</u> , <u>BC052889</u> , <u>NM_016776</u> , <u>NM_016776.1</u> , <u>NM_016776.2</u> , <u>BC004078</u> , <u>BC015070</u> , <u>BC020280</u> , <u>BC060167</u> , <u>NM_016776.3</u>
UniProt ID:	<u>Q7TPV4</u>
Synonyms:	AL024407; AU019902; p67MBP; P160; p160MBP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)