

Product datasheet for **MP208475**

Maz Mouse qPCR Primer Pair (NM_010772)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17188
Forward Sequence:	GGATCACCTCAACAGTCACGTC
Reverse Sequence:	TTCTCCTCGTGTCGGACTGTGT
ACCN:	NM_010772 , NM_010772.1 , BC004758 , BC167175 , BQ572608
Synonyms:	MAZI; Pur-1; PUR1; SAF-1; SAF-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)