

Product datasheet for **MP208494**

Mecom Mouse qPCR Primer Pair (NM_021442)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14013
Forward Sequence:	CCTCCACATTCAAGAGCCATGC
Reverse Sequence:	CTCCAGGCATAGTTGACTCTCG
ACCN:	BC139762 , NM_021442 , NM_021442.2 , BC076620 , BC148544 , BC153133 , NM_021442.3
UniProt ID:	P14404
Synonyms:	D630039M04Rik; Evi-1; Evi1; Jbo; Mds; Mds1; Mds1-Evi1; Prdm3; Znfpr1b1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)