

## Product datasheet for **MP208536**

### Myo1b Mouse qPCR Primer Pair (NM\_010863)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17912
Forward Sequence:	GCCTATGAACCTTGCCTGGAAAG
Reverse Sequence:	GGAGTGCTCTCCACAGGAATC
ACCN:	<a href="#">BC046300</a> , <a href="#">BC054786</a> , <a href="#">NM_010863</a> , <a href="#">NM_010863.2</a> , <a href="#">NM_010863.3</a> , <a href="#">NM_010863.4</a> , <a href="#">BB768105</a> , <a href="#">BC026826</a> , <a href="#">BE950809</a> , <a href="#">BY759345</a>
UniProt ID:	<a href="#">P46735</a>
Synonyms:	AA406997; MIH-L; MMIa
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)