

## Product datasheet for **MP208590**

### Mtmr3 Mouse qPCR Primer Pair (NM\_028860)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74302
Forward Sequence:	AGCAGACTCAGGAACGGACATG
Reverse Sequence:	TCAGATTCCGCACATGGCACAC
ACCN:	<a href="#">BC032166</a> , <a href="#">BC081544</a> , <a href="#">NM_028860</a> , <a href="#">NM_028860.1</a> , <a href="#">NM_028860.2</a> , <a href="#">BC026586</a>
Synonyms:	1700092A20Rik; AI255150; AW557713; FYVE-DSP1; mKIAA0371; ZFYVE10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)