

Product datasheet for **MP208617**

Gnptab Mouse qPCR Primer Pair (NM_001004164)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	432486
Forward Sequence:	GGATTCACATCCAGGAAAGTCCC
Reverse Sequence:	TGTCCTCAGAGTGACGCACCTT
ACCN:	<u>BC060638</u> , <u>NM_001004164</u> , <u>NM_001004164.1</u> , <u>NM_001004164.2</u> , <u>BC060638.1</u> , <u>NM_001004164.3</u>
UniProt ID:	<u>Q69ZN6</u>
Synonyms:	EG432486; mKIAA1208
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)