

Product datasheet for **MP208664**

Ndufb10 Mouse qPCR Primer Pair (NM_026684)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68342
Forward Sequence:	CCATCACCTACTTGACGAAGGC
Reverse Sequence:	TGATGTCTGGCACTCGACGGTA
ACCN:	<u>BC031664</u> , <u>BC043013</u> , <u>NM_026684</u> , <u>NM_026684.1</u> , <u>NM_026684.2</u>
UniProt ID:	<u>Q9DCS9</u>
Synonyms:	22kDa; 0610011B04Rik; PDSW
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)