

Product datasheet for **MP208689**

Noxo1 Mouse qPCR Primer Pair (NM_027988)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71893
Forward Sequence:	TCAGCAGGTAGCCTGGTTTCCA
Reverse Sequence:	CACGGATAGCTCATCAGAGCGA
ACCN:	BC019525 , NM_027988 , NM_027988.1 , NM_027988.2 , NM_027988.3 , NM_027988.4 , BC119527 , BC119528 , BC127235 , BC127236
UniProt ID:	Q8VCM2
Synonyms:	2310034C04Rik; hslt; P41NOX; P41NOXA; P41NOXB; P41NOXC; Snx28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)