

Product datasheet for **MP208713**

Nckap1 Mouse qPCR Primer Pair (NM_016965)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50884
Forward Sequence:	CTCCATTGCCTTTCTCTGCTC
Reverse Sequence:	TGCACAGGAAAGGCTACGATC
ACCN:	NM_016965 , NM_016965.1 , NM_016965.2 , NM_016965.3 , BC003962 , BC032206 , BC054383 , BC060197 , BC127067 , BC148377 , BP764873 , BU055892
UniProt ID:	P28660
Synonyms:	C79304; H19; Hem-2; Hem2; mh19; mKIAA0587; Nap1; p125Nap1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)