

## Product datasheet for **MP208722**

### Necap2 Mouse qPCR Primer Pair (NM\_025383)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66147
Forward Sequence:	GTGTGACAGACTCCAGCAGGTA
Reverse Sequence:	GGTCCTGTAATGCCACGTTGAAG
ACCN:	<u><a href="#">BC037069</a></u> , <u><a href="#">NM_025383</a></u> , <u><a href="#">NM_025383.1</a></u> , <u><a href="#">NM_025383.2</a></u> , <u><a href="#">NM_025383.3</a></u> , <u><a href="#">BB851186</a></u> , <u><a href="#">BK000657</a></u>
UniProt ID:	<u><a href="#">Q9D1J1</a></u>
Synonyms:	1110005F07Rik; AA409105; C78898; NECAP-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)