

## Product datasheet for **MP208843**

### Nrcam Mouse qPCR Primer Pair (NM\_176930)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	319504
Forward Sequence:	GCAGAGTGAACATGACCACACC
Reverse Sequence:	GTGTAGGTTCCGCATCGTCAT
ACCN:	<a href="#">NM_176930</a> , <a href="#">NM_176930.1</a> , <a href="#">NM_176930.2</a> , <a href="#">NM_176930.3</a> , <a href="#">NM_176930.4</a> , <a href="#">BB641673</a> , <a href="#">BC055053</a> , <a href="#">BC060532</a> , <a href="#">BC156235</a> , <a href="#">BC157007</a>
UniProt ID:	<a href="#">Q810U4</a>
Synonyms:	Bravo; C030017F07Rik; C130076O07Rik; mKIAA0343
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)