

Product datasheet for **MP208984**

NdrG1 Mouse qPCR Primer Pair (NM_008681)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17988
Forward Sequence:	TTGCTGTCTGCCATGTGGATGC
Reverse Sequence:	GACTCCAGGAAGCATTTCAGCC
ACCN:	<u>BC015282</u> , <u>BC071235</u> , <u>NM_008681</u> , <u>NM_008681.2</u> , <u>BC141246</u> , <u>BC171968</u> , <u>BY232155</u>
UniProt ID:	<u>Q62433</u>
Synonyms:	CAP43; CMT4D; DRG1; HMSNL; Ndr1; NdrI; NMSL; PROXY1; RTP; TDD5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)