

Product datasheet for **MP209001**

Nipa2 Mouse qPCR Primer Pair (NM_023647)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	93790
Forward Sequence:	CATTCATGCTCCAAAGGAAGAGG
Reverse Sequence:	CATGTCGAGGTCCCACCACAAA
ACCN:	<u>BC038499</u> , <u>NM_023647</u> , <u>NM_023647.1</u> , <u>NM_023647.2</u> , <u>NM_023647.3</u> , <u>NM_023647.4</u> , <u>NM_023647.5</u> , <u>NM_023647.6</u> , <u>BC011415</u> , <u>BC031109</u> , <u>BK001121</u> , <u>BY144304</u>
UniProt ID:	<u>Q9JJC8</u>
Synonyms:	2600017P10Rik; 3830408P04Rik; AB041581
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)