

Product datasheet for **MP209056**

Nfatc2 Mouse qPCR Primer Pair (NM_010899)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18019
Forward Sequence:	ACTTCACAGCGGAGTCCAAGGT
Reverse Sequence:	GGATGTGCTTGTTCCGATACTCG
ACCN:	<u>BC064006</u> , <u>NM_010899</u> , <u>NM_010899.1</u> , <u>NM_010899.2</u> , <u>NM_010899.3</u> , <u>BC064006.1</u> , <u>BP757417</u>
UniProt ID:	<u>Q60591</u>
Synonyms:	AI607462; NF-ATc2; NF-ATp; NFAT1; NFAT1-D; Nfatp
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)