

Product datasheet for **MP209063**

Nfkbib Mouse qPCR Primer Pair (NM_010908)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18036
Forward Sequence:	TCAGCATGAGCCCTTCTGGAT
Reverse Sequence:	CAAGGATGGCTGCTAGATGCAG
ACCN:	<u>BC021938</u> , <u>NM_010908</u> , <u>NM_010908.1</u> , <u>NM_010908.2</u> , <u>NM_010908.3</u> , <u>NM_010908.4</u> , <u>NM_010908.5</u> , <u>BM899229</u>
UniProt ID:	<u>Q60778</u>
Synonyms:	I(Kappa)B(beta); I-kappa-B-beta; Iκ; IKapp; IKappaBbeta; IκB; IκB-B; IκB-beta; IκBβ; NF-kappa-BIB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)