

Product datasheet for **MP209090**

Ncoa2 Mouse qPCR Primer Pair (NM_008678)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17978
Forward Sequence:	TGTCTGATGGCACTCTCGTTGC
Reverse Sequence:	CCATCGCTTGTCAGTCAGATC
ACCN:	<u>BC053387</u> , <u>NM_008678</u> , <u>NM_008678.1</u> , <u>NM_008678.2</u> , <u>NM_008678.3</u> , <u>BP761474</u> , <u>BU703942</u>
UniProt ID:	<u>Q61026</u>
Synonyms:	bHLHe75; D1Ert433e; GRIP-1; Grip1; KAT13C; NCoA-2; SRC-2; TIF2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)