

Product datasheet for **MP209121**

Nr4a3 Mouse qPCR Primer Pair (NM_015743)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18124
Forward Sequence:	ACGCCGAAACCGATGTCAGTAC
Reverse Sequence:	CTCCTGTTGTAGTGGCTCTTTG
ACCN:	<u>NM_015743</u> , <u>NM_015743.1</u> , <u>NM_015743.2</u> , <u>NM_015743.3</u> , <u>BC038520</u> , <u>BC068150</u> , <u>BC128307</u> , <u>BC128308</u> , <u>BC148436</u> , <u>BC156736</u>
UniProt ID:	<u>Q9QZB6</u>
Synonyms:	AI573420; CHN; CSMF; M; MINOR; N; NO; NOR-1; Nor1; TEC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)