

## Product datasheet for **MP209132**

### **Nfxl1 Mouse qPCR Primer Pair (NM\_133921)**

#### **Product data:**

|                   |                                                                                                                                                                                                                                                                                                                                                   |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                                 |
| Gene ID:          | 100978                                                                                                                                                                                                                                                                                                                                            |
| Forward Sequence: | AGCGAATGCCTTGTCACAAGCC                                                                                                                                                                                                                                                                                                                            |
| Reverse Sequence: | CACAAGGTGGACAGTTTCCAGG                                                                                                                                                                                                                                                                                                                            |
| ACCN:             | <u>NM_133921</u> , <u>NM_133921.1</u> , <u>NM_133921.2</u> , <u>BC022652</u> , <u>BC032981</u> , <u>NM_133921.3</u>                                                                                                                                                                                                                               |
| Synonyms:         | 1700012H24Rik; AW538212; D430033A06Rik; DNABF; GCF; Gm1058; TCF9                                                                                                                                                                                                                                                                                  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                       |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage:          | Store at -20 °C.                                                                                                                                                                                                                                                                                                                                  |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                      |



[View online »](#)