

Product datasheet for **MP209137**

Nxt2 Mouse qPCR Primer Pair (NM_172782)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	237082
Forward Sequence:	TGGACAAAGCCACTCTAATCTGG
Reverse Sequence:	CGAGCACTGTAGTTTGGCACTG
ACCN:	<u>BC068166</u> , <u>BC094570</u> , <u>BC118503</u> , <u>NM_172782</u> , <u>NM_172782.1</u> , <u>NM_172782.2</u> , <u>NM_172782.3</u> , <u>NM_172782.4</u> , <u>BB038253</u> , <u>BB851047</u> , <u>BB857403</u> , <u>BC064727</u>
UniProt ID:	<u>Q3UNA4</u>
Synonyms:	6330587F24Rik; AI787442; P15-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)