

Product datasheet for **MP209157**

Nol8 Mouse qPCR Primer Pair (NM_001081350)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70930
Forward Sequence:	ACGGAAGAGACATCCACTCAGG
Reverse Sequence:	TGAACCTGGTGCTGCTCCTT
ACCN:	<u>BC125575</u> , <u>BC125579</u> , <u>NM_001081350</u> , <u>NM_001081350.1</u> , <u>BC026402</u> , <u>BC049166</u> , <u>BC058699</u> , <u>BC094941</u>
UniProt ID:	<u>Q3UHX0</u>
Synonyms:	4921532D18Rik; 5730412B09Rik; D13ErtD548e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)