

## Product datasheet for **MP209212**

### Nudt16l1 Mouse qPCR Primer Pair (NM\_025839)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Gene ID:          | 66911                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | CAACTGCATGCTGTGGAGATCAG                                                                                                                                                                                                                                                                                                                       |
| Reverse Sequence: | ATTGGCAGTGCTGACGAAGGC                                                                                                                                                                                                                                                                                                                         |
| ACCN:             | <u><a href="#">BC025569</a></u> , <u><a href="#">BC028494</a></u> , <u><a href="#">NM_025839</a></u> , <u><a href="#">NM_025839.1</a></u> , <u><a href="#">NM_025839.2</a></u> , <u><a href="#">NM_025839.3</a></u> , <u><a href="#">NM_025839.4</a></u> , <u><a href="#">BQ030979</a></u>                                                    |
| UniProt ID:       | <u><a href="#">Q8VHN8</a></u>                                                                                                                                                                                                                                                                                                                 |
| Synonyms:         | 1110001K21Rik; 5330437I08Rik; Sdos                                                                                                                                                                                                                                                                                                            |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |



[View online »](#)