

Product datasheet for **MP209233**

Ola1 Mouse qPCR Primer Pair (NM_025942)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67059
Forward Sequence:	CTGTGCCAGATGAGAGGTTCTGA
Reverse Sequence:	CCAAGCCCTGCCATTGTGAG
ACCN:	<u>BC011318</u> , <u>BC029207</u> , <u>NM_025942</u> , <u>NM_025942.1</u> , <u>NM_025942.2</u> , <u>BG069761</u> , <u>BX521012</u>
UniProt ID:	<u>Q9CZ30</u>
Synonyms:	2510025G09Rik; 2810405J23Rik; 2810409H07Rik; Gtpbp9; PTD004
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)