

Product datasheet for **MP209235**

Obsl1 Mouse qPCR Primer Pair (NM_178884)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	98733
Forward Sequence:	CAGCTTGAGGATTCTGGCGAGT
Reverse Sequence:	CAAACCTCACGGCGTGTAAGGTC
ACCN:	<u>BC051643</u> , <u>NM_178884</u> , <u>NM_178884.1</u> , <u>NM_178884.2</u> , <u>NM_178884.3</u> , <u>NM_178884.4</u> , <u>NM_178884.5</u> , <u>BC051643.1</u> , <u>BB620298</u> , <u>BC025848</u> , <u>BC056205</u> , <u>BE988113</u> , <u>BI133287</u> , <u>BQ033268</u> , <u>BQ770422</u> , <u>BX512756</u> , <u>BY730281</u>
UniProt ID:	<u>D3YYU8</u>
Synonyms:	AW822216
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)