

Product datasheet for **MP209257**

Olfml2b Mouse qPCR Primer Pair (NM_177068)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	320078
Forward Sequence:	AAGCCAAGGTCTCTGAGGAGGA
Reverse Sequence:	GGCACTTGTCAGTAGGTCTTCC
ACCN:	<u>BC025654</u> , <u>NM_177068</u> , <u>NM_177068.1</u> , <u>NM_177068.2</u> , <u>NM_177068.3</u> , <u>NM_177068.4</u> , <u>BC023743</u>
UniProt ID:	<u>Q3V1G4</u>
Synonyms:	1110018N05Rik; 4832415H08Rik; AI467542
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)