

Product datasheet for **MP209705**

Olfr156 Mouse qPCR Primer Pair (NM_019474)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29846
Forward Sequence:	CTTCGTGCTCATCCTGCTGATG
Reverse Sequence:	AGCAGATGTCCAGGAAGGAGAG
ACCN:	<u>BC101957</u> , <u>BC101961</u> , <u>BC103551</u> , <u>BC104727</u> , <u>NM_019474</u> , <u>NM_019474.1</u> , <u>NM_019474.2</u>
Synonyms:	mOR37b; MOR262-6; Olfr37b; OR37B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)