

Product datasheet for **MP209749**

Olfr205 Mouse qPCR Primer Pair (NM_001011736)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	257881
Forward Sequence:	CCTCCTGGTAGCAATGGCTTATG
Reverse Sequence:	GAACACAGTGCAGAGTCTCTTGG
ACCN:	<u>NM_001011736</u> , <u>NM_001011736.1</u> , <u>BC150839</u>
Synonyms:	GA_x5J8B7U0478-6-340; GA_x54KRFPKG5P-55543875-55542958; MOR182-11P
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)