

Product datasheet for **MP209766**

Olfr23 Mouse qPCR Primer Pair (NM_010970)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18321
Forward Sequence:	CCATGTACCTCACCCTGTCCT
Reverse Sequence:	GAGGTCAGTGAAGGACAAGTTGC
ACCN:	<u>NM_010970</u> , <u>NM_010970.1</u> , <u>BC140105</u> , <u>BC140192</u> , <u>BC141712</u> , <u>BC143029</u> , <u>BC143030</u> , <u>BC146410</u> , <u>BC153106</u>
Synonyms:	MOR135-27; MTPCR50
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)