

Product datasheet for **MP209825**

Olfr320 Mouse qPCR Primer Pair (NM_207230)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216783
Forward Sequence:	CCATCTCTCCTTCCACTGGTGT
Reverse Sequence:	ACTCTGTCCCTTTCCTGAACCC
ACCN:	<u>NM_207230</u> , <u>NM_207230.1</u> , <u>BC104264</u> , <u>BC111096</u> , <u>BC152316</u>
Synonyms:	GA_x6K02T2NKPP-730312-729392; MOR285-4; MOR285-5; Olfr321-ps1; RP23-128D9.8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)