

Product datasheet for **MP209874**

Olfr384 Mouse qPCR Primer Pair (NM_207224)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	193053
Forward Sequence:	GGGAACACTTGATACTGTGGTGC
Reverse Sequence:	GTCCCATAGAACAGTGAGACCAC
ACCN:	<u>NM_207224</u> , <u>NM_207224.1</u> , <u>BC141570</u>
Synonyms:	GA_x5J8B7U2A8F-1-327; GA_x6K02T2P1NL-3773152-3774090; MOR135-5; Olfr384-ps1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)