

Product datasheet for **MP209894**

Olfr417 Mouse qPCR Primer Pair (NM_207137)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	258238
Forward Sequence:	ACCTGGCTCTTTGTGGTCTTCC
Reverse Sequence:	GTGTAGCAGGTCTCAGAGAAGG
ACCN:	<u>NM_207137</u> , <u>NM_207137.1</u> , <u>NM_207137.2</u> , <u>NM_207137.3</u> , <u>BC141574</u>
Synonyms:	GA_x5J8B7W3B3M-312879-313274; GA_x6K02T2P20D-20787051-20786119; MOR267-11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)