

## Product datasheet for **MP210224**

### **Olfr809 Mouse qPCR Primer Pair (NM\_146324)**

#### **Product data:**

|                   |                                                                                                                                                                                                                                                                                                                                                   |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                                 |
| Gene ID:          | 258321                                                                                                                                                                                                                                                                                                                                            |
| Forward Sequence: | CCTTCTCAACCTGCTCCTCACA                                                                                                                                                                                                                                                                                                                            |
| Reverse Sequence: | CAGAGGTGTTGAGCATAGCCAC                                                                                                                                                                                                                                                                                                                            |
| ACCN:             | <u>BC120625</u> , <u>NM_146324</u> , <u>NM_146324.1</u> , <u>BC137816</u>                                                                                                                                                                                                                                                                         |
| Synonyms:         | MOR108-4                                                                                                                                                                                                                                                                                                                                          |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                       |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage:          | Store at -20 °C.                                                                                                                                                                                                                                                                                                                                  |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                      |



[View online »](#)