

Product datasheet for **MP210316**

Olfr924 Mouse qPCR Primer Pair (NM_207560)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	404322
Forward Sequence:	GACTCTGTCTTGCTCAAGCACC
Reverse Sequence:	GAACGTCCTTCATTGGAGCGGA
ACCN:	<u>NM_207560</u> , <u>NM_207560.1</u> , <u>BC141590</u> , <u>BC143111</u> , <u>BC143129</u> , <u>BC143156</u> , <u>BC143208</u>
Synonyms:	MOR171-27P; MOR171-47; Olfr1520-ps1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

