

Product datasheet for **MP210373**

Olfr1274 Mouse qPCR Primer Pair (NM_146263)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	258330
Forward Sequence:	CACCTGAGAACTCAGAGTGCTG
Reverse Sequence:	AAGAGTGGCTGTCTCGTCTCAGA
ACCN:	<u>NM_146263</u> , <u>NM_146263.1</u> , <u>NM_146263.2</u> , <u>BC145937</u>
Synonyms:	GA_x5J8B7W2BV0-3116-4045; MOR228-4; MOR228-5; Olfr1274
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)