

Product datasheet for **MP210380**

Nabp1 Mouse qPCR Primer Pair (NM_028696)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	109019
Forward Sequence:	GAACAGAAGGATAAACTGAGCACC
Reverse Sequence:	GTAGCTGTGGACTGATAGGTCC
ACCN:	BC095967 , BC131923 , BC132532 , NM_028696 , NM_028696.1 , NM_028696.2 , NM_028696.3 , BB609822 , BC027371 , BY089463
UniProt ID:	Q8BGW5
Synonyms:	4930434H03Rik; 4930442A21Rik; 4930488J04Rik; 4933440J18Rik; 5830411E10Rik; AI852561
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)