

Product datasheet for **MP210386**

Oit1 Mouse qPCR Primer Pair (NM_146050)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18300
Forward Sequence:	CATGAGTCCAGTGCGAAACAACG
Reverse Sequence:	CCGTAAGGAAGTTCAGCAGGAG
ACCN:	BC024504 , BC024532 , BC024542 , BC034193 , NM_146050 , NM_146050.1 , NM_146050.2 , BC024636
UniProt ID:	P97805
Synonyms:	2310076N21Rik; AV067083; EF-7; Fam3d
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)