

Product datasheet for **MP210390**

Onecut1 Mouse qPCR Primer Pair (NM_008262)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15379
Forward Sequence:	GAGGATGTGGAAGTGGCTGCAG
Reverse Sequence:	TTGCCTCTGTCCTCCCGTGTT
ACCN:	<u>BC023444</u> , <u>BC024053</u> , <u>NM_008262</u> , <u>NM_008262.1</u> , <u>NM_008262.2</u> , <u>NM_008262.3</u>
UniProt ID:	<u>O08755</u>
Synonyms:	D9Ert423e; Hfh12; HNF-6; Hnf6; Hnf6b; OC-1; Oc1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)