

Product datasheet for **MP210409**

Tmem259 Mouse qPCR Primer Pair (NM_001003949)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216157
Forward Sequence:	AGCCGCCTATTGCTGGATGAGT
Reverse Sequence:	GCCATCCACATGCTGACAAAGC
ACCN:	<u>NM_001003949</u> , <u>NM_001003949.1</u> , <u>NM_001003949.2</u> , <u>NM_001003949.3</u> , <u>BB645384</u> , <u>BC005494</u> , <u>BC052787</u> , <u>BC056654</u> , <u>BC080302</u> , <u>BC100297</u> , <u>BC156275</u> , <u>BC157059</u> , <u>NM_001003949.4</u>
UniProt ID:	<u>Q8CIV2</u>
Synonyms:	AC087114.7gm2; mFLJ00277; ORF61; R32184_3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)