

Product datasheet for **MP210419**

Oprl1 Mouse qPCR Primer Pair (NM_011012)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18389
Forward Sequence:	GCTCAGCACAAGTGGAGGATGA
Reverse Sequence:	GGCTGTAGCAGACAGAGATGATC
ACCN:	BC050885 , BC051982 , NM_011012 , NM_011012.1 , NM_011012.2 , NM_011012.3 , NM_011012.4 , NM_011012.5 , BU708032 , BY267034
UniProt ID:	P35377
Synonyms:	KOR-3; KOR3; LC13; LC132; mo; MOR; MOR-C; morc; NOP; O; Oprl; ORGC; ORL1; XO; XOR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)