

Product datasheet for **MP210440**

Orc4 Mouse qPCR Primer Pair (NM_011958)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26428
Forward Sequence:	AGTCCACATAGCAACCTCTTTGG
Reverse Sequence:	GTGGTCTTCTGATCCTCGTG
ACCN:	BC015072 , BC019748 , NM_011958 , NM_011958.1 , NM_011958.2 , NM_011958.3 , BB620704 , BB672099 , BM117783 , NM_011958.4
UniProt ID:	O88708
Synonyms:	mMmORC4; Orc4l; Orc4P
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)