

Product datasheet for **MP210456**

Otx1 Mouse qPCR Primer Pair (NM_011023)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18423
Forward Sequence:	GTGGCACTCAAGATCAACCTGC
Reverse Sequence:	GGAGACGACTTCTTCTTGACCG
ACCN:	BC057105 , BC058354 , NM_011023 , NM_011023.1 , NM_011023.2 , NM_011023.3 , BB652255
UniProt ID:	P80205
Synonyms:	A730044F23Rik; j; jv
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

