

Product datasheet for **MP210460**

Osgepl1 Mouse qPCR Primer Pair (NM_028091)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72085
Forward Sequence:	TTACGTGCTGGCTTGGCGTTT
Reverse Sequence:	CGCGGTACTTTTATGGCAGCTTC
ACCN:	BC058172 , NM_028091 , NM_028091.1 , NM_028091.2 , NM_028091.3 , BC030671 , BC038910 , BC051211 , BY780319
UniProt ID:	Q6PEB4
Synonyms:	2610001M19Rik; AA416452
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)