

Product datasheet for **MP210500**

Oxsr1 Mouse qPCR Primer Pair (NM_133985)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	108737
Forward Sequence:	GTTCTGGACGAGCCAACCATTG
Reverse Sequence:	GCAATCTGCACAGAGCCGTCTT
ACCN:	BC060645 , NM_133985 , NM_133985.1 , NM_133985.2 , BC020039 , BC036148 , BC046453 , BC057992
UniProt ID:	Q6P9R2
Synonyms:	2210022N24Rik; 2810422B09Rik; AI462649; AW209236; mKIAA1101; Osr1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)