

Product datasheet for **MP210503**

Oxgr1 Mouse qPCR Primer Pair (NM_001001490)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	239283
Forward Sequence:	TGTATCTGACCAGCCTCCCGTT
Reverse Sequence:	CAGGTGAGGAAGAGAATGCTGC
ACCN:	<u>BC119196</u> , <u>BC119200</u> , <u>NM_001001490</u> , <u>NM_001001490.1</u> , <u>NM_001001490.2</u> , <u>NM_001001490.3</u>
UniProt ID:	<u>Q6IYF8</u>
Synonyms:	Cysltr3; Gm299; Gpr80; Gpr99; P2Y15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)