

Product datasheet for **MP210510**

Osbp1a Mouse qPCR Primer Pair (NM_207530)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64291
Forward Sequence:	CCTGTCAGCATCCTTAGTGAGG
Reverse Sequence:	CAAAGGAGTGTGCTGTCACAGC
ACCN:	<u>NM_207530</u> , <u>NM_207530.1</u> , <u>NM_207530.2</u> , <u>NM_207530.3</u> , <u>BB667982</u> , <u>BC015075</u> , <u>BC031735</u> , <u>BC046466</u> , <u>BC057194</u> , <u>BC076637</u> , <u>BC092247</u> , <u>BC148511</u> , <u>BC153110</u> , <u>BG092905</u> , <u>BY248497</u> , <u>BY249673</u> , <u>BY715608</u>
UniProt ID:	<u>Q91XL9</u>
Synonyms:	G430090F17Rik; Gm753; ORP-1; Osbp1b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)