

## Product datasheet for **MP210522**

### Pak2 Mouse qPCR Primer Pair (NM\_177326)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	224105
Forward Sequence:	GCTGTAGTGACAGAGGAAGATGA
Reverse Sequence:	TCACCAACTGGAGCAGGAATGG
ACCN:	<u><a href="#">BC086650</a></u> , <u><a href="#">NM_177326</a></u> , <u><a href="#">NM_177326.1</a></u> , <u><a href="#">NM_177326.2</a></u> , <u><a href="#">NM_177326.3</a></u> , <u><a href="#">BC029195</a></u> , <u><a href="#">BC043712</a></u> , <u><a href="#">BC063068</a></u>
UniProt ID:	<u><a href="#">Q8CIN4</a></u>
Synonyms:	5330420P17Rik; A130002K10Rik; AI836325; D16ErtD269e; mKIAA4182; PAK-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)