

Product datasheet for **MP210526**

Pak7 Mouse qPCR Primer Pair (NM_172858)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	241656
Forward Sequence:	CTCATCCTCAAGGACACTCGTAC
Reverse Sequence:	GGAGGACTGAAGACCATCTGTG
ACCN:	<u>NM_172858, NM_172858.1, NM_172858.2, BC099417, BC148654, BC156990</u>
UniProt ID:	<u>Q8C015</u>
Synonyms:	2900083L08Rik; 6430627N20; PAK-5; PAK-7; Pak5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)